Lovibond® Water Testing

Tintometer® Group

Product Care

This product has been produced to the highest standard.
Prior to use please store in a cool place at around 10 - 15°C.

Large variations in storage temperature will cause some condensation inside the tube & eventual dehydration of the agar.

Samples from water & fluids

Remove the paddle & cap assembly by pulling gently while twisting clockwise, take care not to touch the agar surface.

Immerse the paddle in the sample to be tested or expose to a spray or running fluid and ensure the agar is covered.

After allowing excess liquid to drain replace the paddle & cap in the tube and ensure this is pushed fully home.

The dipslide should now immediately be incubated at 30 - 35°C, check the results after 1 - 2 days for bacteria and up to 5 days for yeast & moulds.

Dipslides

GB

Surface samples

Remove the paddle & cap assembly by pulling gently while twisting clockwise, take care not to touch the agar surface.

Press each agar side of the paddle in different locations on the desired surfaces, replace the paddle & cap in the tube and ensure this pushed fully home.

The dipslide should now immediately be incubated at 30 - 35°C, check the results after 1 - 2 days for bacteria and up to 5 days for yeast & moulds.

Constant temperature incubation is important to ensure results are as accurate as possible, if the temperature is lower than recommended microbial growth may be slow, leading to false interpretation. If you are unsure it is advisable to continue incubation for a further period.

Interpretation

When you are ready the growth may be compared with the charts on the rear of this leaflet.

Water & fluids

If the sample contains very high levels of contamination this may lead to confluent growth due to saturation etc. Please compare with a new sterile slide if uncertain.

Surfaces

The area of each agar is 10 cm². It has been shown the 50 % of viable organisms are recovered from surfaces tested. To increase the sensitivity sweep the surface area during testing.

Please dispose of used dipslides carefully,

they should be autoclaved, incinerated or soaked in fresh 5 - 10 % bleach solution overnight.

Code	Type	Description	Growth
56B010110	D001 TTC/TTC	TVC (Aerobic bacteria)	Red
56B010210	D002 TTC/MALT	TVC/Fungi & Moulds	Red/Green, Black, White: Beige (Yeast)
56B010310	D003 TTC/ROSE	TVC/Fungi & Moulds	Red/Green, Black, White: Beige (Yeast)
56B010410	D004 TTC/MAC	TVC/Coliforms	Red/Red or Yellow
56B010510	D005 TTC/E.COLI	TVC/ Coliforms & <i>E.coli</i>	Red/ Red or Yellow (Coliforms) or Blue (<i>E.coli</i>) or Cream (<i>Paeruginosa</i>)
56B010610	D006 MAC/PDM	Coliforms/ <i>Pseudomonas</i>	Red/Green (<i>Paeruginosa</i>) or Cream (<i>P.spp</i>)
56B010710	D007 TTC/PDM	TVC/ <i>Paeruginosa</i>	Red/ Green (<i>Paeruginosa</i>) or Cream (<i>P.spp</i>)
56B010810	D008 SRB Tube Test	Sulphate Reducing Bacteria	Black Diffusion
56B010910	D009 NRB Tube Test	Nitrite Reducing Bacteria	Pink Diffusion
56B011110	D011 R2A/R2A	Aerobic bacteria	Red

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are registered trademarks of the
The Tintometer Limited
Tintometer Inc.