

# Lovibond® Water Testing

## Tintometer® Group



### Step 1

Please open the instruction „HexLoad\_de\_gb.pdf“ in the „HexLoad“ dir.  
Print the instruction and have it at hand when updating.

### Step 2

To update your SpectroDirect, MultiDirect, PoolDirect, PC9in1, TurbiDirect, PCPECTRO II or MaxiDirect open the program HEXLoad.

### Step 3

Additional every Product has its own firmware. You will find it in the folder „Firmware“ fitting to your product.

### Step 4

Update your device as described in the instruction „HexLoad\_de\_gb.pdf“.

### Step 5

Download the current version of the instruction for your product from our website.

## Instruction Software-Update

# Lovibond® Water Testing

## Tintometer® Group



## Changelog

### New regarding MD600/MultiDirect:

- New Methods Bromine PP,  $\text{ClO}_2$  PP,  $\text{H}_2\text{O}_2$  LR and  $\text{H}_2\text{O}_2$  HR
- New way to handle M265 Nitrate tube test

---

### New regarding PM620/PoolDirect:

- New Method  $\text{H}_2\text{O}_2$  HR