

The ZIP file need to be unpacked and the folder "Exchange\_Method\_Profile" needs to be copied to the main directory of USB Stick. From this USB stick you can load the methods onto your XD7000 and XD7500.

To import methods from USB stick please follow these steps:

1. Select "General Setup" from main menu.
2. Select "Exchange methods / profiles".
3. Select "Import from USB memory device"
4. All methods in the USB will be pre-selected for import, unselect the methods which are not required.
5. Press "F4" to import the selected methods.
6. The imported methods will appear at the end of methods list in Concentration menu.

**Note:** if a method number of imported methods already exists, the instrument will ask to overwrite the existing method or to cancel the import.

Methods descriptions are also included as pdf. in this ZIP folder.

For information about reagents, please visit our website.