

Safety data sheet according to 1907/2006/EC, Article 31

Printing date 10.11.2023

Version number 27 (replaces version 26)

Revision: 10.11.2023

SECTION 1: Identification of the substance/mixture and of the company/undertaking

· **1.1 Product identifier**

· **Product name:** Hydrogenperoxide LR

· **Catalog number:** 00512381, 512380BT, 512381BT, 00512389BT, 4512380BT, 4512381BT

· **1.2 Relevant identified uses of the substance or mixture and uses advised against**

· **Application of the substance / the preparation:** Reagent for water analysis

· **1.3 Details of the supplier of the safety data sheet**

· **Supplier:**

Tintometer GmbH
Schleefstraße 8-12
44287 Dortmund
Made in Germany
www.lovibond.com

phone: +49 (0)231 94510-0
e-mail: sales@lovibond.com

The Tintometer Limited
Lovibond® House
Sun Rise Way
Amesbury
Wiltshire SP4 7GR
United Kingdom

phone : +44 1980 664800
e-mail: SDS@lovibond.uk

· **Informing department:**

e-mail: sds@lovibond.com
Product Safety Department

· **1.4 Emergency telephone number:**

+44 1235 239670
Languages: English

SECTION 2: Hazards identification

· **2.1 Classification of the substance or mixture**

· **Classification according to Regulation (EC) No 1272/2008**



GHS08 health hazard

STOT RE 1 H372 Causes damage to the thyroid through prolonged or repeated exposure. Route of exposure: Oral.

· **2.2 Label elements**

· **Labelling according to Regulation (EC) No 1272/2008**

The product is classified and labelled according to the GB CLP regulation.

· **Hazard pictograms**



GHS08

· **Signal word** Danger

· **Hazard-determining components of labelling:**

potassium iodide

· **Hazard statements**

H372 Causes damage to the thyroid through prolonged or repeated exposure. Route of exposure: Oral.

· **Precautionary statements**

P264 Wash hands thoroughly after handling.

(Contd. on page 2)

— GB —

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 10.11.2023

Version number 27 (replaces version 26)

Revision: 10.11.2023

Product name: Hydrogenperoxide LR

(Contd. of page 1)

P314 Get medical advice/attention if you feel unwell.

2.3 Other hazards

The main intake pathways of potassium iodide are: inhalation of dust and solution aerosols, as well as oral ingestion.

Results of PBT and vPvB assessment

This mixture does not contain any substances that are assessed to be persistent, bioaccumulative and toxic (PBT) or very persistent and very bioaccumulative (vPvB), according to the criteria given in Annex XIII of Regulation (EC) No. 1907/2006.

Determination of endocrine-disrupting properties

The product does not contain substances with endocrine disrupting properties.

SECTION 3: Composition/information on ingredients

3.2 Mixtures
Description: Mixture of organic and inorganic compounds

Dangerous components:

CAS: 7681-11-0	potassium iodide	 STOT RE 1, H372	20–30%
EINECS: 231-659-4			
Reg.nr.: 01-2119966161-40-XXXX			

Additional information For the wording of the listed hazard phrases refer to section 16.

SECTION 4: First aid measures

4.1 Description of first aid measures
General information Instantly remove any clothing soiled by the product.

After inhalation Supply fresh air; consult doctor in case of symptoms.

After skin contact

Instantly wash with water and soap and rinse thoroughly.

If skin irritation continues, consult a doctor.

After eye contact

Rinse opened eye for several minutes under running water (at least 15 min). If symptoms persist, consult doctor.

After swallowing

Rinse out mouth and then drink 1-2 glasses of water.

In case of persistent symptoms consult doctor.

4.2 Most important symptoms and effects, both acute and delayed:

irritating effects possible

after swallowing and inhalation:

absorption

general feeling of sickness

after absorption of large amounts:

headache

thirst

sickness

vomiting

gastric pain

diarrhoea

drop in blood pressure

cardiovascular disorders

weakness

4.3 Indication of any immediate medical attention and special treatment needed:

Absorption: in case of iodine hypersensitivity, even after relatively low doses, acute respiratory and cardiovascular disorders (possibly shock), skin and mucous membrane reactions possible. (GESTIS)

Symptoms of poisoning may even occur after several hours.

SECTION 5: Firefighting measures

5.1 Extinguishing media
Suitable extinguishing agents Use fire fighting measures that suit the environment.

5.2 Special hazards arising from the substance or mixture

The product is not combustible.

Formation of toxic gases is possible during heating or in case of fire.

Can be released in case of fire:

(Contd. on page 3)

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 10.11.2023

Version number 27 (replaces version 26)

Revision: 10.11.2023

Product name: Hydrogenperoxide LR

(Contd. of page 2)

nitrous gases
 Nitrogen oxides (NO_x)
 Sulphur oxides (SO_x)
 Phosphorus oxides (P_xO_x)
 Dipotassium oxide
 Sodium oxide
 Hydrogen iodide (HI)
 smoke of metal oxide

· **5.3 Advice for firefighters**

· **Protective equipment:**

Wear self-contained breathing apparatus.
 Wear full protective suit.

· **Additional information**

Collect contaminated fire fighting water separately. It must not enter drains.
 Dispose of fire debris and contaminated fire fighting water in accordance with official regulations.
 Ambient fire may liberate hazardous vapours.

SECTION 6: Accidental release measures

· **6.1 Personal precautions, protective equipment and emergency procedures**

· **Advice for non-emergency personnel:**

Wear protective equipment. Keep unprotected persons away.
 Ensure adequate ventilation

· **Advice for emergency responders:** Protective equipment: see section 8

· **6.2 Environmental precautions:** Do not allow product to reach sewage system or water bodies.

· **6.3 Methods and material for containment and cleaning up:**

Ensure adequate ventilation.
 Collect mechanically.
 Dispose of contaminated material as waste according to item 13.

· **6.4 Reference to other sections**

See Section 8 for information on personal protection equipment.
 See Section 13 for information on disposal.

SECTION 7: Handling and storage

· **7.1 Precautions for safe handling**

· **Advice on safe handling:** Provide suction extractors if dust is formed.

· **Hygiene measures:**

Take off immediately all contaminated clothing.
 Wash hands during breaks and at the end of the work.
 Do not eat, drink or smoke when using this product.

· **7.2 Conditions for safe storage, including any incompatibilities**

· **Requirements to be met by storerooms and containers:** Store in cool location.

· **Information about storage in one common storage facility:** Store away from oxidising agents.

· **Further information about storage conditions:**

Store in a locked cabinet or with access restricted to technical experts or their assistants.
 Protect from heat and direct sunlight.

Store in cool, dry conditions in well sealed containers.

Protect from the effects of light.

Protect from humidity and keep away from water.

This product is hygroscopic.

· **Recommended storage temperature:** 20°C +/- 5°C

· **7.3 Specific end use(s)** No further relevant information available.

GB
 (Contd. on page 4)

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 10.11.2023

Version number 27 (replaces version 26)

Revision: 10.11.2023

Product name: Hydrogenperoxide LR

(Contd. of page 3)

SECTION 8: Exposure controls/personal protection

· 8.1 Control parameters

· Components with limit values that require monitoring at the workplace:

CAS: 7631-95-0 sodium molybdate

WEL (Great Britain)	Short-term value: 10 mg/m ³ Long-term value: 5 mg/m ³ as Mo
---------------------	---

· Regulatory information WEL (Great Britain): EH40/2020

· DNELs

CAS: 7681-11-0 potassium iodide

Oral	DNEL	0.01 mg/kg /bw/d (Consumer / acute / systemic effects)
		0.01 mg/kg /bw/d (Consumer / long-term / systemic effects)
Dermal	DNEL	1 mg/kg /bw/d (Worker / long-term /systemic effects)
		1 mg/kg /bw/d (Consumer / long-term / systemic effects)
Inhalative	DNEL	0.07 mg/m ³ (Worker / long-term /systemic effects)
		0.035 mg/m ³ (Consumer / long-term / systemic effects)

· PNECs

CAS: 7681-11-0 potassium iodide

PNEC	0.007 mg/l (Fresh water)
PNEC	0.075 mg/kg (Aquatic intermittent release)
	0.007 mg/kg /sediment (Fresh water sediment)

· **Additional information:** The lists that were valid during the compilation were used as basis.

· 8.2 Exposure controls

· Engineering measures:

Technical measures and appropriate working operations should be given priority over the use of personal protective equipment. See item 7.

· Individual protection measures, such as personal protective equipment

Protective clothing should be selected specifically for the workplace, depending on concentration and quantity of the hazardous substances handled.

· Eye/face protection

Safety glasses

use against the effects of fumes / dust

Use safety glasses that have been tested and approved in accordance with government standards such as EN 166.

· Hand protection

Preventive skin protection by use of skin-protecting agents is recommended.

After use of gloves apply skin-cleaning agents and skin cosmetics.

· Material of gloves

nitrile rubber, NBR

Recommended thickness of the material: ≥ 0.11 mm

· Penetration time of glove material

Value for the permeation: Level = 1 (< 10 min)

The exact break trough time has to be found out by the manufacturer of the protective gloves and has to be observed.

· Other skin protection (body protection):

Protective work clothing.

· Breathing equipment:

Use breathing protection against the effects of fumes/dust/aerosol.

· Recommended filter device for short term use:

Filter P3

· Environmental exposure controls

Do not allow product to reach sewage system or water bodies.

SECTION 9: Physical and chemical properties

· 9.1 Information on basic physical and chemical properties

· Physical state	Solid.
· Form:	Tablets
· Colour:	Light grey
· Odour:	Odourless
· Odour threshold:	Not applicable.

(Contd. on page 5)

GB

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 10.11.2023

Version number 27 (replaces version 26)

Revision: 10.11.2023

Product name: Hydrogenperoxide LR

(Contd. of page 4)

· Melting point/Freezing point:	Not determined.
· Boiling point or initial boiling point and boiling range	Not determined.
· Flammability	The product is not combustible.
· Explosive properties:	Product is not explosive.
· Lower and upper explosion limit	
Lower:	Not applicable.
Upper:	Not applicable.
· Flash point:	Not applicable.
· Auto-ignition temperature:	Not applicable (solid).
· Decomposition temperature:	Not determined.
· pH (12.5 g/l) at 20°C	5.5
· Kinematic viscosity	Not applicable (solid).
· Solubility	
· Water:	Soluble
· Partition coefficient n-octanol/water (log value)	Not applicable (mixture).
· Vapour pressure:	Not applicable (solid).
· Density and/or relative density	
· Density:	Not determined.
· Relative density:	Not determined.
· Relative gas density	Not applicable (solid).
· Particle characteristics	Not determined.
· 9.2 Other information	
· Information with regard to physical hazard classes	
· Corrosive to metals	Void
· Other safety characteristics	
· Oxidising properties:	none
· Additional information	
· Solids content:	100.0 %

SECTION 10: Stability and reactivity

- **10.1 Reactivity** No further relevant information available.
- **10.2 Chemical stability** Stable at ambient temperature (room temperature).
- **10.3 Possibility of hazardous reactions**
 Reacts with acids, alkalis and oxidizing agents
 --> forms heat
 Reacts with peroxides
 Reacts with halogenated compounds
- **10.4 Conditions to avoid** Strong heating (decomposition)
- **10.5 Incompatible materials:** No further relevant information available.
- **10.6 Hazardous decomposition products:** see section 5

SECTION 11: Toxicological information

- **11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008**
- **Acute toxicity** Based on available data, the classification criteria are not met.

· **LD/LC50 values that are relevant for classification:**

CAS: 7681-11-0 potassium iodide

Oral	LD50	2779 mg/kg (rat)
Dermal	LD50	3160 mg/kg (rabbit)
	NOAEL	0.01 mg/kg /bw/d (human) organ: Thyroid

- **Skin corrosion/irritation** Based on available data, the classification criteria are not met.
- **Serious eye damage/irritation** Based on available data, the classification criteria are not met.
- **Information on components:** CAS 6283-63-2: DPD may cause allergic skin reaction
- **Respiratory or skin sensitisation** Based on available data, the classification criteria are not met.

(Contd. on page 6)

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 10.11.2023

Version number 27 (replaces version 26)

Revision: 10.11.2023

Product name: Hydrogenperoxide LR

(Contd. of page 5)

Information on components:

CAS 6283-63-2: Sensitization possible in predisposed persons.

The following applies to iodides in general: Sensitization possible at predisposed persons.

- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.

- **Carcinogenicity** Based on available data, the classification criteria are not met.

- **Reproductive toxicity** Based on available data, the classification criteria are not met.

Information on components:

OECD 414: Teratogenicity testing

OECD 473: Mutagenicity testing

OECD 471, 474, 476, 487: Germ cell mutagenicity testing

CAS: 7681-11-0 potassium iodide

OECD 471 (negative) (Bacterial Reverse Mutation Test - Ames test)

OECD 476 (negative) (In Vitro Mammalian Cell Gene Mutation Test)

Mouse (lymphoma L5178Y cells)

- **STOT (specific target organ toxicity) -single exposure** Based on available data, the classification criteria are not met.

- **STOT (specific target organ toxicity) -repeated exposure**

Causes damage to the thyroid through prolonged or repeated exposure. Route of exposure: Oral.

- **Aspiration hazard** Based on available data, the classification criteria are not met.

Information on likely routes of exposure

"Main routes of exposure:

At workplaces, intake of potassium iodide (KI) is most likely to occur via the respiratory tract.

Outside the workplace, iodides are ingested with food (essential) and sometimes with medications.

Respiratory tract: KI can be inhaled as dust or aerosol from solutions. Inhalation studies were conducted with particulate aerosols containing sodium iodide using various animal species (monkey, mouse, sheep). Rapid and effective absorption via the respiratory tract was observed. This is also assumed for KI as its solubility is comparable.

Skin: From tests on volunteers who had an aqueous KI solution applied to their forearms (12.5 cm²), the amount of iodine absorbed was estimated at 0.1%. Absorption through the skin is therefore considered to be of little relevance.

Gastrointestinal tract: Soluble iodide is absorbed almost entirely via the gastrointestinal tract. This has been proven by results of studies with KI on adult volunteers." [GESTIS]

Additional toxicological information:

In case of an acute molybdenum(VI) intoxication: diarrhoea, anaemia, fatigue, loss of appetite. Toxic effect on liver and kidneys after high doses.

CAS: 7681-11-0 potassium iodide

(source: GESTIS)

Main Toxic Effects:

Acute: Irritation to the eyes, skin and airways, disturbance of thyroid function, cardiovascular effects, metabolic disturbances.

Chronic: Disturbance of thyroid function, systemically conditioned skin damage and inflammation of the mucous membranes.

Further Information (GESTIS, Merck):

Small amounts of iodine are essential for the body. However, long-term overdoses of iodine lead to disturbances in the thyroid function (hypo- and/or hyperthyroidism, possibly accompanied by thyroiditis). The effects are very complex.

Furthermore, symptoms of chronic iodine poisoning (iodine toxicosis, "iodism") can occur following intake of high doses of predisposed persons. They mainly consist of systemically conditioned irritation/inflammatory changes to the mucous membranes and skin.

Iodide crosses the placenta and, when administered (orally) to pregnant women in very high doses, can lead to hypothyroidism and/or goiter in the fetus with deaths from tracheal compression

11.2 Information on other hazards

- **Endocrine disrupting properties** The product does not contain substances with endocrine disrupting properties.

Other information

Other dangerous properties can not be excluded.

According to the information available to us, the chemical, physical and toxicological properties of the substances mentioned in Chapter 3 have not been thoroughly investigated.

GB
(Contd. on page 7)

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 10.11.2023

Version number 27 (replaces version 26)

Revision: 10.11.2023

Product name: Hydrogenperoxide LR

(Contd. of page 6)

SECTION 12: Ecological information

· 12.1 Toxicity

· Aquatic toxicity:

CAS: 7681-11-0 potassium iodide

EC50 7.5 mg/l/48h (Daphnia magna) (OECD 202)
MerckLC50 3780 mg/l/96h (rainbow trout) (OECD 203)
Merck

· Other information:

Toxic for fish:

molybdenum compounds in general: > 25 mg/l

· **12.2 Persistence and degradability** No further relevant information available.· **12.3 Bioaccumulative potential** No further relevant information available.· **12.4 Mobility in soil** No further relevant information available.

· 12.5 Results of PBT and vPvB assessment

This mixture does not contain any substances that are assessed to be persistent, bioaccumulative and toxic (PBT) or very persistent and very bioaccumulative (vPvB), according to the criteria given in Annex XIII of Regulation (EC) No. 1907/2006.

· **12.6 Endocrine disrupting properties** The product does not contain substances with endocrine disrupting properties.

· 12.7 Other adverse effects

Depending on the concentration, phosphorus and/or nitrogen compounds may contribute to the eutrophication of water supplies. Avoid transfer into the environment.

· Water hazard:

Do not allow product to reach ground water, water bodies or sewage system, even in small quantities. Danger to drinking water if even extremely small quantities leak into soil.

SECTION 13: Disposal considerations

· 13.1 Waste treatment methods

· **Recommendation** Must not be disposed of together with household garbage. Do not allow product to reach sewage system.

· European waste catalogue

16 05 06* laboratory chemicals, consisting of or containing hazardous substances, including mixtures of laboratory chemicals

· Uncleaned packagings:

· **Recommendation:** Disposal must be made according to official regulations.· **Recommended cleaning agent:** Water, if necessary with cleaning agent.

SECTION 14: Transport information

· 14.1 UN number or ID number

· ADR, IMDG, IATA Void

· 14.2 UN proper shipping name

· ADR, IMDG, IATA Void

· 14.3 Transport hazard class(es)

· ADR, IMDG, IATA
· Class Void

· 14.4 Packing group

· ADR, IMDG, IATA Void

· **14.5 Environmental hazards:** Not applicable.· **14.6 Special precautions for user** Not applicable.

· 14.7 Maritime transport in bulk according to IMO instruments

Not applicable.

· Transport/Additional information:

Not dangerous according to the above specifications.

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 10.11.2023

Version number 27 (replaces version 26)

Revision: 10.11.2023

Product name: Hydrogenperoxide LR

(Contd. of page 7)

SECTION 15: Regulatory information

· 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

· Poisons Act UK

· Regulated explosives precursors

None of the ingredients is listed.

· Regulated poisons

None of the ingredients is listed.

· Reportable explosives precursors

None of the ingredients is listed.

· Reportable poisons

None of the ingredients is listed.

· Regulation (EU) 2019/1148 on the marketing and use of explosives precursors not regulated

· Regulation (EU) No 649/2012 concerning the export and import of hazardous chemicals (PIC)

None of the ingredients is listed.

· Regulation (EC) No 1334/2000 setting up a Community regime for the control of exports of dual-use items and technology:

None of the ingredients is listed.

· Regulation (EC) No 273/2004 on drug precursors

None of the ingredients is listed.

· Regulation (EC) No 111/2005 laying down rules for the monitoring of trade between the Community and third countries in drug precursors

None of the ingredients is listed.

· Regulation (EC) No 1005/2009 on substances that deplete the ozone layer:

None of the ingredients is listed.

· REGULATION (EU) 2019/1021 on persistent organic pollutants (POP)

None of the ingredients is listed.

· LIST OF SUBSTANCES SUBJECT TO AUTHORISATION (ANNEX XIV)

None of the ingredients is listed.

· Substances of very high concern (SVHC) according to REACH, Article 57

This product does not contain any substances of very high concern above the legal concentration limit of $\geq 0.1\%$ (w / w).

· Substances of very high concern (SVHC) according to UK REACH

This product does not contain any substances of very high concern above the legal concentration limit of $\geq 0.1\%$ (w / w).

· Directive 2012/18/EU (SEVESO III):

· Named dangerous substances - ANNEX I None of the ingredients is listed.

· Information about limitation of use: Employment restrictions concerning young persons must be observed (94/33/EC).

· 15.2 Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

SECTION 16: Other information

These data are based on our present knowledge. However, they shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

This Safety Data Sheet is in compliance with Regulation (EC) No 1907/2006, Article 31 as amended by Regulation (EU) 2020/878.

· Training hints Provide adequate information, instruction and training for operators.

· Relevant phrases

H372 Causes damage to organs through prolonged or repeated exposure.

· Abbreviations and acronyms:

OECD: Organisation for Economic Co-operation and Development

STOT: specific target organ toxicity

SE: single exposure

RE: repeated exposure

EC50: half maximal effective concentration

IC50: half maximal inhibitory concentration

NOEL or NOEC: No Observed Effect Level or Concentration

(Contd. on page 9)

GB

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 10.11.2023

Version number 27 (replaces version 26)

Revision: 10.11.2023

Product name: Hydrogenperoxide LR

(Contd. of page 8)

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)

RID: Règlement international concernant le transport des marchandises dangereuses par chemin de fer (Regulations Concerning the International Transport of Dangerous Goods by Rail)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

DNEL: Derived No-Effect Level (UK REACH)

PNEC: Predicted No-Effect Concentration (UK REACH)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

SVHC: Substances of Very High Concern

vPvB: very Persistent and very Bioaccumulative

STOT RE 1: Specific target organ toxicity (repeated exposure) – Category 1

Sources

Data arise from safety data sheets, reference works and literature.

ECHA: European Chemicals Agency <http://echa.europa.eu>

GESTIS- Stoffdatenbank (Substance Database, Germany)

* Data compared to the previous version altered.

GB