

Lovibond® Water Testing

Tintometer® Group



Step 1

Please open the instruction „HexLoad_de_gb.pdf“ in the „HexLoad“ dir.
Print the instruction and have it at hand when updating.

Step 2

To update your SpectroDirect, MultiDirect, PoolDirect, PC9in1, TurbiDirect, PCPECTRO II or MaxiDirect open the program HEXLoad.

Step 3

Additional every Product has its own firmware. You will find it in the folder „Firmware“ fitting to your product.

Step 4

Update your device as described in the instruction „HexLoad_de_gb.pdf“.

Step 5

Download the current version of the instruction for your product from our website.

Instruction Software-Update