# Lovibond® Colour Measurement

# Tintometer® Group



### IP Units Filter Set



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

### IP Units Filter Set

A set of four glass conformance filters for checking the PFXi instruments in IP (method B) units

Part Number: 137660

#### **Technical Data**

 Glass Standard Values
 Approximate values

 - 0.75
 - 1.25

 - 2.0
 - 3.5

Filters

## **Delivery Scope**

- IP 0.75 filter
- IP 1.25 filter
- ── IP 2.0 filter
- —— IP 3.5 filter
- Certificate of conformity