

# Lovibond® Water Testing Tintometer® Group



## DI 20 Incubator

Compliant with WHO guidelines



Our new DI 20 incubator has been designed to allow the user to reach global drinking water standards for testing, even in difficult environments.

The DI 20 incubator is a high quality, portable, lightweight incubator for microbiological testing of water using plates. It is the only incubator of its type to come with heating and cooling as standard ensuring the samples are incubated at the correct temperature no matter what ambient conditions.

Its innovative design also includes a small petri dish holder utilising our re-usable stainless steel micro plates, the dishes allow each incubator to perform up to 20 tests.

Its sleek design and intuitive pictographic buttons ensure that it can be used with little or no training and its LED lights allow users to see at a glance how far through a cycle the incubation is, and that the temperature is stable.

- Compliant with WHO guidelines " Guidelines for drinking water quality 4th edition."
- Compact design & lightweight
- Heating & Cooling fitted as standard
- Intuitive, unique pictographic buttons

Part Number: 56B000714

### Industry

NGO

### Technical Data

<b>Interfaces</b>	USB
<b>Power Supply</b>	12 V
<b>Temperature Range</b>	20 °C - 47 °C
<b>Temperature Tolerance</b>	± 0.5 °C
<b>Environmental Conditions</b>	5 °C - 50 °C
<b>Compliance</b>	CE
<b>Dimensions</b>	116 x 165 x 116 mm

### Delivery Scope

- DI 20
- Power leads
- 50 Petri dishes
- Power Supply

## Accessories

Title	Part Number
Power cable set	136300
USB Cable 2.0 A - Micro B plug	190630
Petri dish holder	197139
Power supply, 12 V / 48 W	19803550
Battery cable	19803555
Petri dishes, 10 pc.	400855

---

### Tintometer GmbH

Lovibond® Water Testing  
Schleefstraße 8-12  
44287 Dortmund  
Tel.: +49 (0)231/94510-0  
sales@lovibond.com  
www.lovibond.com  
Germany

### The Tintometer Limited

Lovibond House  
Sun Rise Way  
Amesbury, SP4 7GR  
Tel.: +44 (0)1980 664800  
Fax: +44 (0)1980 625412  
sales@lovibond.uk  
www.lovibond.com  
UK

### Tintometer China

9F, SOHO II C.  
No.9 Guanghualu,  
Chaoyang District,  
Beijing, 100020  
Customer Care China Tel.: 4009021628  
Tel.: +86 10 85251111 Ext. 330  
Fax: +86 10 85251001  
chinaoffice@tintometer.com  
www.lovibond.com  
China

### Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,  
Lebuh Nilam 2, Bandar Bukit Tinggi,  
Klang, 41200, Selangor D.E  
Tel.: +60 (0)3 3325 2285/6  
Fax: +60 (0)3 3325 2287  
lovibond.asia@tintometer.com  
www.lovibond.com  
Malaysia

### Tintometer Brazil

Caixa Postal: 271  
CEP: 13201-970  
Jundiaí – SP  
Tel.: +55 (11) 3230-6410  
sales@lovibond.us  
www.lovibond.com.br  
Brazil

### Tintometer Inc.

6456 Parkland Drive  
Sarasota, FL 34243  
Tel: 941.756.6410  
Fax: 941.727.9654  
sales@lovibond.us  
www.lovibond.us  
USA

### Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2<sup>nd</sup>, 3<sup>rd</sup> & 4<sup>th</sup> Floor  
Sanathnagar Industrial Estate,  
Hyderabad, 500018  
Telangana  
Tel: +91 (0) 40 23883300  
Toll Free: 1 800 599 3891/ 3892  
indiaoffice@lovibond.in  
www.lovibondwater.in  
India

### Tintometer Spain

Postbox: 24047  
08080 Barcelona  
Tel.: +34 661 606 770  
sales@tintometer.es  
www.lovibond.com  
Spain