

Lovibond® Colour Measurement Tintometer® Group



Filtre individuel (certifié)

Code: 109970



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Filtre individuel (certifié)

This part number stands for a single filter supplied in other filter sets. Please note that each individual filter has its own part number. Please contact sales@tintometer.com for the correct number for your filter.

Fiche technique

Glass Standard Values	Préciser l'instrument, l'échelle et la valeur nominale à la commande.
Filters	1

Contenu de livraison

- 1 filter usually supplied in another set
- Certificate of conformity