

# Lovibond<sup>®</sup> Colour Measurement Tintometer<sup>®</sup> Group



## Padrão líquido Platinum-cobalto número de cor 15 - 500 ml

Artigo No: 134160



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

### Padrão líquido Platinum-cobalto número de cor 15 - 500 ml

Garrafa de 500 ml, certificada segundo o sistema QS ISO 9001 com um prazo de validade de 12 meses.

Rastreabilidade completa de acordo com normas reconhecidas internacionalmente

Cada garrafa é fornecida com certificados completos, incluindo ficha de dados de segurança.