

# Lovibond<sup>®</sup> Colour Measurement Tintometer<sup>®</sup> Group



## Saybolt 色标 +25

货号: 134070



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

### Saybolt 色标 +25

500 ml 瓶子，经过 ISO/IEC 17025 认证，保质期为 12 个月。

完全可追溯性符合国际公认标准

每个瓶子均随附完整证书，包括安全数据表。