

Lovibond® Water Testing Tintometer® Group



Sewage / Effluent Test Kit

Test all important parameters for sewage/effluent analysis



- Easy to use
- Futuristic design
- RAPID tablets fast dissolving
- Highest accuracy

Part Number: 56K701220

Industry

Industries Others

Application

Cooling Water | Disinfection Control

Accessories

Title	Part Number
CHECKIT® Disc pH (Universal)	146130
Plastic shaker tube, graduated 100 ml, with lid	385130A
Titration jar with cap, plastic, 60 mL	56A006701
Turbidity tube (two piece) - 5-500 JTU (NTU)	56A010801
Turbidity tube (one piece) - 30-400 NTU	56A607101

Delivery Scope

- Plastic case
- Checkit Comparator
- 4 - 10 pH Universal Checkit Disc
- Universal pH tablets
- Permanganat tablets
- Acidifying Sewage Effluent tablets
- 2 plastic vials 10 mL with lid
- Shaker tube with lid 100 mL
- Syringe plastic, 20 mL
- Titration jar with cap, plastic, 60 mL
- Turbidity tube
- Instruction manual