

Lovibond® Water Testing Tintometer® Group



Chlorine Dioxide HR Test Strips

For fast determination with chemical reaction



- Hazard free
- Rapid
- Always to hand
- Semi quantitative results

Part Number: 56S000550

Industry

Industries Others | Municipalities

Application

Disinfection Control

Chlorine Dioxide HR Test Strips

Test strips offer quick, clean results without the fuss of reagents or tablets. Their dip and compare method makes them a popular choice for many applications where speed and simplicity are most important.

Measuring Range

Test Name	Measuring Range
Chlorine Dioxide LR	0.5 - 100 mg/L ClO ₂