

Lovibond® Water Testing Tintometer® Group



Hydrogen Peroxide Strips

For fast determination with chemical reaction



- Hazard free
- Rapid
- Always to hand
- Semi quantitative results

Part Number: 56S000850

Industry

Industries Others

Application

Cooling Water | Disinfection Control

Hydrogen Peroxide Strips

Test strips offer quick, clean results without the fuss of reagents or tablets. Their dip and compare method makes them a popular choice for many applications where speed and simplicity are most important.

Measuring Range

Test Name	Measuring Range
Hydrogen Peroxide test strips	0.5 - 100 mg/L H ₂ O ₂