

# Lovibond® Colour Measurement Tintometer® Group



## Filtro específico para el usuario

Referencia No: 13109981



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

### Filtro específico para el usuario

A single filter of a type not usually produced as part of a standard filter set. Please contact [sales@tintometer.com](mailto:sales@tintometer.com) with the details of your requirement to check that it is technically possible to create the required filter.

### Datos técnicos

Glass Standard Values Según corresponda

Filters 1

## Volumen de suministro

- User specified glass filter
- Certificate of conformity