


# Lovibond® Colour Measurement Tintometer® Group



## Juego de filtros EP rojo

Referencia No: 139410



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

### Juego de filtros EP rojo

A set of two glass conformance filters for checking PFXi instruments in European Pharmacopoeia Red scale.

### Datos técnicos

<b>Glass Standard Values</b>	- R2
	- R6

---

<b>Filters</b>	2
----------------	---

### Volumen de suministro

- EP Red 2 filter
- EP Red 6 filter
- Certificate of Conformity