

Lovibond® Colour Measurement Tintometer® Group



Lot de filtres Ph.Eur marron

Code: 139430



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Lot de filtres Ph.Eur marron

A set of two glass conformance filters for checking PFXi instruments in European Pharmacopoeia Brown scale.

Fiche technique

Glass Standard Values	B2, B5
Filters	2

Contenu de livraison

- EP Brown 2 filter
- EP Brown 5 filter
- Certificate of conformity