

# Lovibond® Colour Measurement Tintometer® Group



## Lot de filtres Ph.Eur marron/jaune

Code: 139440



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

### Lot de filtres Ph.Eur marron/jaune

A set of two glass conformance filters for checking PFXi instruments in European Pharmacopoeia Brown/Yellow scale.

### Fiche technique

<b>Glass Standard Values</b>	BY3, BY5
<b>Filters</b>	2

### Contenu de livraison

- EP Brown/Yellow 3 filter
- EP Brown/Yellow 5 filter
- Certificate of Conformity