

Lovibond® Colour Measurement Tintometer® Group



Bêta-carotène

Code: 139520



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Bêta-carotène

A glass conformance filter to check the PFXi range of instruments accuracy in measuring beta Carotene.

Fiche technique

Glass Standard Values	45 ppm avec un trajet optique de 10 mm
Filters	1

Contenu de livraison

- 45ppm filter
- Certificate of Conformity