

Lovibond® Colour Measurement Tintometer® Group



Lot de filtres Lovibond® rouge, jaune, neutre (PFXi 995/950/880)

Code: 139590



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Lot de filtres Lovibond® rouge, jaune, neutre (PFXi 995/950/880)

A set of ten glass conformance filters for checking PFXi instruments (PFXi 995/950/880) in Lovibond RYBN scale.

Fiche technique

Glass Standard Values	Valeurs approximatives :
	- 0.1R 0.5Y 0.0N
	- 0.5R 1.3Y 0.1N
	- 1.0R 2.6Y 0.1N
	- 1.3R 3.3Y 0.5N
	- 2.0R 5.7Y 0.1N
	- 2.5R 8.6Y 0.1N
	- 3.2R 15.0Y 0.1N
	- 4.0R 22.0Y 0.1N
	- 5.0R 34.0Y 0.1N
	- 7.1R 51.0Y 0.5N

Contenu de livraison

- Lovibond RYBN 0.1R 0.5Y 0.0N filter
- Lovibond RYBN 0.5R 1.3Y 0.1N filter
- Lovibond RYBN 1.0R 2.6Y 0.1N filter
- Lovibond RYBN 1.3R 3.3Y 0.5N filter
- Lovibond RYBN 2.0R 5.7Y 0.1N filter
- Lovibond RYBN 2.5R 8.6Y 0.1N filter
- Lovibond RYBN 3.2R 15.0Y 0.1N filter
- Lovibond RYBN 4.0R 22.0Y 0.1N filter
- Lovibond RYBN 5.0R 34.0Y 0.1N filter
- Lovibond RYBN 7.1R 51.0Y 0.5N filter
- Certificate of Conformity

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Allemagne

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
Royaume-Uni

Tintometer China

Room 1001, China Life Tower
16 Chaoyangmenwai Avenue,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 App. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
Chine

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Malaisie

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
Brésil

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
États-Unis

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
Inde

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
Espagne