

Lovibond® Colour Measurement Tintometer® Group



Filtre personnalisé (Fx)

Code: 169980



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Filtre personnalisé (Fx)

A single filter of a type not usually produced as part of a standard filter set. Please contact sales@tintometer.com with the details of your requirement to check that it is technically possible to create the required filter.

Fiche technique

Glass Standard Values Selon le cas

Filters	1
---------	---

Contenu de livraison

- ☐ User specified glass filter
- ☐ Certificate of conformity