

Lovibond® Colour Measurement Tintometer® Group



Filtro singolo (certificato)

N. ordine: 109970



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Filtro singolo (certificato)

This part number stands for a single filter supplied in other filter sets. Please note that each individual filter has its own part number. Please contact sales@tintometer.com for the correct number for your filter.

Dati tecnici

Glass Standard Values	specify instrument, scale and nominal value when ordering.
Filters	1

Fornitura

- 1 filter usually supplied in another set
- Certificate of conformity