

# Lovibond<sup>®</sup> Colour Measurement Tintometer<sup>®</sup> Group



## Standard liquido numero di colore Platino-Cobalto 5 - 500 ml

N. ordine: 134140



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

### Standard liquido numero di colore Platino-Cobalto 5 - 500 ml

Bottiglia da 500 ml, certificato secondo il sistema QS ISO 9001 con una durata di conservazione di 12 mesi.

Tracciabilità completa secondo standard riconosciuti a livello internazionale

Ogni bottiglia viene fornita con certificati completi, compresa la scheda tecnica di sicurezza.