

# Lovibond® Colour Measurement Tintometer® Group



## Numero di colore AOCS-Tintometer 0.4R 2.0Y (5¼")

N. ordine: 134240



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

## Numero di colore AOCS-Tintometer 0.4R 2.0Y (5¼")

Bottiglia da 500 ml, certificato secondo il sistema QS ISO 9001 con una durata di conservazione di 12 mesi.

Tracciabilità completa secondo standard riconosciuti a livello internazionale

Ogni bottiglia viene fornita con certificati completi, compresa la scheda tecnica di sicurezza.