

Lovibond[®] Colour Measurement Tintometer[®] Group



Numero di colore ASTM <0,5

N. ordine: 134290



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

Numero di colore ASTM <0,5

Bottiglia da 500 ml certificata ISO/IEC 17025 con una durata di conservazione di 12 mesi.

Tracciabilità completa secondo standard riconosciuti a livello internazionale

Ogni bottiglia viene fornita con certificati completi, compresa la scheda tecnica di sicurezza.