

# Lovibond<sup>®</sup> Colour Measurement Tintometer<sup>®</sup> Group



## ASTM-kleurcode <0,5

Bestelnr.: 134290



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

### ASTM-kleurcode <0,5

500-ml-fles, gecertificeerd volgens ISO/IEC 17025 met een houdbaarheid van 12 maanden.

Volledige traceerbaarheid volgens internationaal erkende normen

Elke fles wordt geleverd met volledige certificaten, inclusief veiligheidsinformatieblad.