

Lovibond® Colour Measurement Tintometer® Group



Filterset EP rood

Bestelnr.: 139410



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Filterset EP rood

A set of two glass conformance filters for checking PFXi instruments in European Pharmacopoeia Red scale.

Technische gegevens

| | |
|-----------------------|------|
| Glass Standard Values | - R2 |
| | - R6 |

| | |
|---------|---|
| Filters | 2 |
|---------|---|

Leveromvang

- EP Red 2 filter
- EP Red 6 filter
- Certificate of Conformity