

# Lovibond® Colour Measurement Tintometer® Group



## Conjunto de filtros EP Amarelo

Artigo No: 139420



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

### Conjunto de filtros EP Amarelo

A set of two glass conformance filters for checking PFXi instruments in European Pharmacopoeia Red scale.

### Dados Técnicos

Glass Standard Values Y2, Y4

Filters 2

### Texto do Escopo de Entrega

- EP Yellow 2 filter
- EP Yellow 4 filter
- Certificate of Conformity