

Lovibond® Colour Measurement Tintometer® Group



Шкала Lovibond® красный, желтый, нейтральный (Fх)

Номер заказа: 139209



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Шкала Lovibond® красный, желтый, нейтральный (Fх)

A glass conformance filter for checking Model Fx instruments in Lovibond RYBN scale.

Технические данные

Glass Standard Values	Approximate Values
	- 4.5R 19.2Y 0.8N

Filters	1
---------	---

Объем поставки

- Lovibond RYBN 4.5R 19.2Y 0.8N filter
- Certificate of Conformity