

# Lovibond® Colour Measurement Tintometer® Group



## Набор фильтров для шкалы EP коричневый/ желтый

Номер заказа: 139440



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

### Набор фильтров для шкалы EP коричне- вый/желтый

A set of two glass conformance filters for checking PFXi instruments in European Pharmacopoeia Brown/Yellow scale.

#### Технические данные

Glass Standard Values	- BY3 - BY5
-----------------------	----------------

Filters	2
---------	---

### Объем поставки

- EP Brown/Yellow 3 filter
- EP Brown/Yellow 5 filter
- Certificate of Conformity