

# Lovibond® Colour Measurement Tintometer® Group



## Набор фильтров для шкалы EP зеленый/ желтый

Номер заказа: 139450



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

### Набор фильтров для шкалы EP зеленый/жел- тый

A set of two glass conformance filters for checking PFXi instruments in European Pharmacopoeia Green/Yellow scale.

#### Технические данные

Glass Standard Values	- GY2
	- GY5

Filters	2
---------	---

### Объем поставки

- EP Green/Yellow 2 filter
- EP Green/Yellow 5 filter
- Certificate of Conformity