

Lovibond® Colour Measurement Tintometer® Group



AOCS-Tintometer Color 3.0R 28.0Y (5¼")



- Certified according to ISO 17025 / ISO 17034 under UKAS accreditation
- Full Traceability to International Standards.
- Tamper evident security packaging

Part Number: 134280

Note: Reference liquid standards are assigned a nominal value and Lovibond always aims to match this value as closely as possible. The value achieved and the expected performance tolerances are reported on the Certificate supplied with the standard.

AOCS-Tintometer Color 3.0R 28.0Y (5¼")

ASTM (ASTM D6045, D1500) Color Reference Standards are ideal for routine verification of color measurement instruments, verifying test data and ensuring good correlation of results between instruments or laboratories. Each standard is issued with a Certificate of Calibration which confirms measurement values, a printed SDS and is UKAS Certified.

Supplies 500 ml.