

# Lovibond® Colour Measurement Tintometer® Group



## EP Brown Filter Set

Part Number: 139430



Conformance filters are assigned a nominal value and the company always endeavors to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

### EP Brown Filter Set

A set of two glass conformance filters for checking PFXi instruments in European Pharmacopoeia Brown scale.

### Technical Data

<b>Glass Standard Values</b>	- B2 - B5
------------------------------	--------------

<b>Filters</b>	2
----------------	---

### Delivery Scope

- EP Brown 2 filter
- EP Brown 5 filter
- Certificate of conformity