

Lovibond[®] Colour Measurement Tintometer[®] Group



Índice de color Saybolt 0

Referencia No: 134050



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

Índice de color Saybolt 0

Frasco de 500 ml, certificado según la norma ISO/IEC 17025 con una durabilidad de 12 meses.

Trazabilidad absoluta de conformidad con los estándares reconocidos a nivel internacional

Cada frasco se suministra con todos los certificados, incl. hoja de datos de seguridad.