

# Lovibond<sup>®</sup> Colour Measurement Tintometer<sup>®</sup> Group



## Liquide de référence Indice cobalt-platine 10 - 500 ml

Code: 134150



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

### Liquide de référence Indice cobalt-platine 10 - 500 ml

Flacon de 500 ml, certifié conforme au système d'assurance qualité ISO 9001, peut être conservé pendant 12 mois.

Traçabilité complète, conforme aux normes internationales reconnues  
Chaque flacon est fourni avec des certificats complets, y compris fiche technique de sécurité.