

# Lovibond® Colour Measurement Tintometer® Group



## Lot de filtres Ph.Eur rouge

Code: 139410



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

### Lot de filtres Ph.Eur rouge

A set of two glass conformance filters for checking PFXi instruments in European Pharmacopoeia Red scale.

### Fiche technique

Glass Standard Values R2, R6

Filters 2

### Contenu de livraison

- EP Red 2 filter
- EP Red 6 filter
- Certificate of Conformity