

# Lovibond<sup>®</sup> Colour Measurement Tintometer<sup>®</sup> Group



## Número de cor Saybolt 0

Artigo No: 134050



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

### Número de cor Saybolt 0

Frasco de 500 ml certificado segundo ISO/IEC 17025 com um prazo de validade de 12 meses.

Rastreabilidade completa de acordo com normas reconhecidas internacionalmente

Cada garrafa é fornecida com certificados completos, incluindo ficha de dados de segurança.