

Lovibond® Water Testing Tintometer® Group



Pooltester Chlorine HR/pH in box

The fast assistants for Pool Testing



High measuring accuracy

Results are shown with high measurement accuracy. They have a certainty and can counteract if necessary.

Intuitive

Self-explanatory is the operation of the tester. And if you are ever unsure: A look at the packaging and the intuitively understandable instructions are sufficient to eliminate the ambiguity.

Ergo mixed

The ergonomically optimized lid ensures a secure grip.

- Easy to use
- Futuristic design
- RAPID tablets fast dissolving
- Highest accuracy

Part Number: 151601

Industry

Marine Industry | Pool Domestic | Pool Services

Application

Disaster Management | Pool Water Control

Pooltester Chlorine HR/pH in box

The pool tester enables the simultaneous measurement of the water treatment agent and the pH value.

Measuring Range

Test Name	Measuring Range
Chlorine	0.5 - 6.0 mg/L Cl ₂
pH	6.8 - 8.2 pH

Delivery Scope

- Pooltester in a sturdy plastic box
- Tablet reagents for 20 tests
- Instruction manual