

Lovibond® Water Testing

Tintometer® Group



Pooltester Chlorine LR/pH in carton

The fast assistants for Pool Testing



- Easy to use
- RAPID tablets fast dissolving
- Highest accuracy

Part Number: 151613

High measuring accuracy

Results are shown with high measurement accuracy. They have a certainty and can counteract if necessary.

Intuitive

Self-explanatory is the operation of the tester. And if you are ever unsure: A look at the packaging and the intuitively understandable instructions are sufficient to eliminate the ambiguity.

Ergo mixed

The ergonomically optimized lid ensures a secure grip.

Industry

Pool Domestic

Application

Cooling Water | Disaster Management | Disinfection Control | Drinking Water Treatment | Food and Beverage | Pool Water Control

Pooltester Chlorine LR/pH in carton

A step towards a sustainable future: the Pooltester Cl LR, pH is now available in recyclable packaging.

The Pooltester enables the simultaneous measurement of chlorine content and pH value: fill in the pool water, add Rapid tablets and compare the colour change with the scale. This makes it easy for pool owners to keep track of their pool's water maintenance.

Measuring Range

| Test Name | Measuring Range |
|-----------|--------------------------------|
| Chlorine | 0.1 - 3.0 mg/L Cl ₂ |
| pH | 6.8 - 8.2 pH |

Delivery Scope

- Pooltester
- Tablet reagents for 20 tests
- Instruction manual

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Germany

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
UK

Tintometer China

9F, SOHO II C.
No.9 Guanghualu,
Chaoyang District,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 Ext. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
China

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebu Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Malaysia

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
Brazil

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
USA

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
India

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
Spain