

# Lovibond® Water Testing

## Tintometer® Group



### Nessleriser 2250 with stand, Daylight unit and nessler cylinders DB 420

More than 400 different test discs available



- More than 400 different test discs available
- Compensation for coloured and turbid samples
- Guaranteed constancy of the coloured glass standards
- Integrated prism

Part Number: 172040

#### Usable in different applications

The Nessleriser 2250 is versatile: it is used in swimming pools as well as in research centres, universities or general drinking water treatment.

#### Daylight lamp

So that you can always measure safely and stably under unfavourable visibility conditions, this package comes with a mains daylight unit, so the lighting conditions always remain the same no matter where you are measuring.

#### Integrated prism

The integrated prism visually brings the glass standards of the colour disc and the coloured sample into view.

#### Industry

Chemical Industry | Food and Beverage Industry | Industries Others | Municipalities | NGO | Oil and Gas | Pharmaceutical Industry | Power and Energy

### Nessleriser 2250 with stand, Daylight unit and nessler cylinders DB 420

The Nessleriser 2250 is a versatile colorimetric system for colour measurement and water analysis. Comfortable handling, no compromise on accuracy and reproducibility: The Nessleriser 2250 can be used for many applications. Swimming pools, research or drinking water treatment are just a few examples.

#### Accessories

Title	Part Number
Daylight cabinet, mains operated	171020
Stand for Upgrade Kit	172180
Pair Nessler tubes, 113 mm path length, 50 ml, AF 306	353060
Pair Nessler tubes, 113 mm path length with plungers, 50 ml, AF 306/P	353080
AF 315 1 Special Nessler tube (determination of Oxygen with colour disc NOE)	353150
Pair Nessler tubes, 250 mm path length with plungers, DB 420	354200