

# Lovibond® Water Testing Tintometer® Group



## Ozone 3/67/0.1-1.0 mg/l Kit

More than 400 different test discs available



- More than 400 different test discs available
- Compensation for coloured and turbid samples
- Guaranteed constancy of the coloured glass standards
- Integrated prism

Part Number: 56K016401

### Usable in different applications

The Comparator 2000+ is versatile: it is used in swimming pools as well as in research centers, universities or general drinking water treatment.

### Daylight lamp

So that you can always measure safely and stably under unfavorable visibility conditions, we recommend the use of the (optional) mains or battery-powered daylight lamp. So the lighting conditions always remain the same. No matter where you are measuring.

### Various analyses possible

For the Comparator 2000+, more than 400 color discs and numerous reagents are available. Your combination is guaranteed. And if you have to change them once, it works without opening the comparator.

### Integrated prism

The integrated prism visually brings the glass standards of the color disc and the colored sample into view.

### Industry

Chemical Industry | Food and Beverage Industry | Industries Others | Marine Industry | Oil and Gas | Pharmaceutical Industry | Pool Public | Power and Energy

### Application

Drinking Water Treatment

### Ozone 3/67/0.1-1.0 mg/l Kit

The Comparator 2000+ is a versatile colorimetric system for water analysis. Comfortable handling, no compromise on accuracy and reproducibility: The Comparator 2000+ can be used for many applications. Swimming pools, research or drinking water treatment are just a few examples.

### Measuring Range

| Test Name | Measuring Range               | Chemical Method |
|-----------|-------------------------------|-----------------|
| Ozone T   | 0.1 - 1.0 mg/L O <sub>3</sub> | DPD             |