

Lovibond® Water Testing Tintometer® Group



Copper Strips 0.05-2 mg/l (25 pc.)

For fast determination with chemical reaction



- Hazard free
- Rapid
- Always to hand
- Semi quantitative results

Part Number: 56S000625

Industry

Industries Others

Application

Cooling Water

Copper Strips 0.05-2 mg/l (25 pc.)

Test strips offer quick, clean results without the fuss of reagents or tablets. Their dip and compare method makes them a popular choice for many applications where speed and simplicity are most important.

Measuring Range

Test Name

Copper-Test Strips

Measuring Range

0.05 - 2 mg/L Cu