

Lovibond® Water Testing

Tintometer® Group



Disinfection Chlorine Test Kit

For monitoring of chlorinated water systems



- Maintenance before and after disinfection works
- Monitoring the correct chlorine concentration & avoiding damages
- Determination Chlorine, pH & more disinfection parameters
- Complete Kit for on-site testing

Referencia No: 56K701130

Chlorination is the process of adding chlorine to water to prevent the spread of infection and disease (including Legionnaires' Disease).

Chlorinated water systems require careful monitoring to ensure the correct conditions are achieved and maintained during the chlorination process. Failure to achieve the correct chlorination levels can lead to unsuccessful disinfections while elevated levels can lead to potentially corrosive conditions.

This Disinfection Chlorine Test Kit has been designed to measure low and high levels of free chlorine. Test strips allow quick easy checking at Sentinel outlets.

Industria

Otras industrias

Aplicaciones

Agua de refrigeración | Control de desinfección

Disinfection Chlorine Test Kit

For commercial buildings, especially dormant premises, industry organisations recommend in many countries regular maintenance procedures for hot & cold-water services to prevent infection with bacteria, biofilms & material damage through misdosing.

The Lovibond® Disinfection Chlorine Test Kit assists you with testing instruments and reagents for monitoring the pH level, chlorine and other important parameters for safe water supply.

Rango de medición

| Test Name | Rango de medición |
|---------------|-------------------------------|
| Cloro (libre) | 1 - 300 mg/L Cl ₂ |
| Cloro HR | 10 - 200 mg/L Cl ₂ |
| pH | 7 - 14 pH |

Volumen de suministroo

- Plastic case
- pH strips, pH 7-14, plastic
- Chlorine HR strips, 10-200 mg/L
- Free Chlorine Buffer FCL1
- Free Chlorine LR Titrant FCL2
- Free Chlorine HR Titrant FCL3
- Syringe plastic, 20 mL
- Titration jar with cap, plastic, 60 mL
- Instruction manual

Accesorios

| Título | Referencia No |
|--|---------------|
| Jeringuilla, plástico, 20 mL | 56A006501 |
| Recipiente de valoración con tapa, plástico, 60 mL | 56A006701 |

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Alemania

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
Reino Unido

Tintometer China

9F, SOHO II C.
No.9 Guanghualu,
Chaoyang District,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 Ext. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
China

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Malasia

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
Brasil

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
Estados Unidos de América

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
India

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
España

Nos reservamos el derecho de cambios técnicos

Printed in Germany

Lovibond® and Tintometer® are Trademarks of the Tintometer Group of Companies