

# Lovibond® Water Testing

## Tintometer® Group



## Sewage / Effluent Test Kit

Test all important parameters for sewage/effluent analysis



- Easy to use
- Futuristic design
- RAPID tablets fast dissolving
- Highest accuracy

Referencia No: 56K701220

### Industria

Otras industrias

### Aplicaciones

Agua de refrigeración | Control de desinfección

## Volumen de suministro

- Plastic case
- Checkit Comparator
- 4 - 10 pH Universal Checkit Disc
- Universal pH tablets
- Permanganat tablets
- Acidifying Sewage Effluent tablets
- 2 plastic vials 10 mL with lid
- Shaker tube with lid 100 mL
- Syringe plastic, 20 mL
- Titration jar with cap, plastic, 60 mL
- Turbidity tube
- Instruction manual

## Accesorios

Título	Referencia No
Disco CHECKIT® pH (universal)	146130
Recipiente de plástico, graduado 100 ml, con tapa	385130A
Recipiente de valoración con tapa, plástico, 60 mL	56A006701
Análisis de turbidez 5-500 JTU	56A010801
Turbidity tube (one piece) - 30-400 NTU	56A607101

---

### Tintometer GmbH

Lovibond® Water Testing  
Schleefstraße 8-12  
44287 Dortmund  
Tel.: +49 (0)231/94510-0  
sales@lovibond.com  
www.lovibond.com  
Alemania

### The Tintometer Limited

Lovibond House  
Sun Rise Way  
Amesbury, SP4 7GR  
Tel.: +44 (0)1980 664800  
Fax: +44 (0)1980 625412  
sales@lovibond.uk  
www.lovibond.com  
Reino Unido

### Tintometer China

9F, SOHO II C.  
No.9 Guanghualu,  
Chaoyang District,  
Beijing, 100020  
Customer Care China Tel.: 4009021628  
Tel.: +86 10 85251111 Ext. 330  
Fax: +86 10 85251001  
chinaoffice@tintometer.com  
www.lovibond.com  
China

### Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,  
Lebuh Nilam 2, Bandar Bukit Tinggi,  
Klang, 41200, Selangor D.E  
Tel.: +60 (0)3 3325 2285/6  
Fax: +60 (0)3 3325 2287  
lovibond.asia@tintometer.com  
www.lovibond.com  
Malasia

### Tintometer Brazil

Caixa Postal: 271  
CEP: 13201-970  
Jundiaí – SP  
Tel.: +55 (11) 3230-6410  
sales@lovibond.us  
www.lovibond.com.br  
Brasil

### Tintometer Inc.

6456 Parkland Drive  
Sarasota, FL 34243  
Tel: 941.756.6410  
Fax: 941.727.9654  
sales@lovibond.us  
www.lovibond.us  
Estados Unidos de América

### Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2<sup>nd</sup>, 3<sup>rd</sup> & 4<sup>th</sup> Floor  
Sanathnagar Industrial Estate,  
Hyderabad, 500018  
Telangana  
Tel: +91 (0) 40 23883300  
Toll Free: 1 800 599 3891/ 3892  
indiaoffice@lovibond.in  
www.lovibondwater.in  
India

### Tintometer Spain

Postbox: 24047  
08080 Barcelona  
Tel.: +34 661 606 770  
sales@tintometer.es  
www.lovibond.com  
España