

Lovibond® Water Testing

Tintometer® Group



Water Safety Kit Combined Duo

Uitgerust voor alle vereisten



You need safe, fast and reliable test results in difficult situations like disasters or potable water monitoring in rural environments without laboratories? You immediate certainty about the quality of potable water to save health and lives? Take all you need with you and combine microbiological methods and simple chemical methods of analysis. Take the smaller Water Safety Kit Chemical in one hand and the bigger Water Safety Kit Microbiology! on your back in the optional backpack. You are able to test different sorts of bacteria on-site due to the heating and cooling function of our DI 20 petri dish incubator.

- Emergency response kits for disaster situations and water monitoring outback
- Equipment for chemical and microbiological analysis in two cases

Bestelnr.: 56K681252

Industrie

Chemische industrie | Energieleveranciers | Farmaceutische industrie | Gemeenten | Ngo | Olie-industrie | Openbaar zwembad | Overige industrieën | Privézwembad | Scheepvaart | Voedingsmiddelen- en drankenindustrie | Zwembadservice

Applikation

Afvalwaterzuivering | Behandeling drinkwater | Controle desinfecteermiddel | Controle zwembadwater | Ketelwater | Koelwater | Zuivering- vervuild water

Water Safety Kit Combined Duo

Combine both microbiological and physio chemical analysis of some of the most critical parameters to determine the suitability of water for drinking. Take the extended complete equipment including 2 incubators for bacteria tests with you everywhere you need immediate and reliable test result.

With their rugged, waterproof cases the kits are easy to transport. Even non technical personal is able to use it due to pictoral instructions.

User

Non Profit Organizations (NGO's) | Medical research and development | Institutes, Universities, Schools | World Health Organization (WHO)

Meetbereik

Test Name	Meetbereik	Chemische methode
Aërobe bacteriën (TVC)/ colibacteriën		
Aërobe bacteriën (TVC)/ E.Coli, coliformen		
Ammonium T	0.02 - 1 mg/L N	Indofenolblauw
Chlorine T	0.01 - 6.0 mg/L Cl ₂ ^{a)}	DPD
Conductivity	0 - 20 mS/cm	
Nitraat T	0.08 - 1 mg/L N	Zinkreductie / NED
pH	0 - 14	
Turbidity	0.01 - 1100 NTU	

Leveromvang

- 1 case Water Safety Kit Microbiology Duo
- 1 case Water Safety Kit Chemical

Accessoires

Omschrijving	Bestelnr.
Stroomkabelset	136300
Incubator Battery, Lead/Acid	192345
Satz Trübungsstandards T-CAL (<0,1, 20, 200, 800 NTU)	194150
Flesjes	363590
Universal container with spoon and lid, 30 ml	366239
Plastic beker, gegradueerd, 100 ml, met deksel	385130
Plastic beker, gegradueerd, 30 ml, met deksel	385132
Petrischaaltjes, 10 stuks	400855
ValidCheck Chloor 1,5 mg/l	48105510
Roerstaaf en poederlepel	56A006601
Rugzak voor Water Safety Kits	56A01400
Media Pads and Dispenser (Pack of 200)	56A016330
Maatbeker	56A018501
Membraanfiltratie	56A023801
Pincet	56A604201
Membraanfiltratie	56B000710
D004 TTC/MAC Dipslides (Macconkey)	56B010410
D005 TTC/Chromogeen Dipslides	56B010510
KS3-pH 10 bufferoplossing 65 ml	56L000365
KS6-pH 4 bufferoplossing 65 ml	56L000665
KS9-pH 7 bufferoplossing 65 ml	56L000965
KS16-1413US Standaardgeleidbaarheidsoplossing	56L001665
Membrane Lauryl Sulphate Broth, 38.1 gr.	56P069740
Pump (35 ml stroke)	56X610301
Pipette, disposable, 1.5 ml	BI-K10060
Marker Pen	BI-K14087

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Duitsland

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
Verenigd Koninkrijk

Tintometer China

9F, SOHO II C.
No.9 Guanghualu,
Chaoyang District,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 Ext. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
China

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Maleisië

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
Brazilië

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
Verenigde Staten

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
India

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
Spanje