


# Lovibond® Water Testing

## Tintometer® Group



### Water Safety Kit Microbiology Duo

Uitgerust voor alle vereisten



Disaster situation or potable water monitoring in demanding situations: You need certainty and reliability for test results in everywhere. The Water Safety Kit Microbiology Duo contains all the equipment necessary to perform the microbiological analysis component for safe drinking water. Due to the heating and cooling function of the DI 20 incubator you can test different sorts of bacteria - with both incubators simultaneously! Perform filtration with sterilisable equipment including sampler, filter flask & vacuum pump on-site - sterilisation equipment not included. With the optional backpack you can take the kit everywhere with you.

- Testing coliforms and E. coli bacteria according to WHO standards
- Testing different forms of bacteria simultaneously with 2 incubators
- DI 20 with heating and cooling as standard

Bestelnr.: 56K681255

#### Industrie

Ngo | Overige industrieën

#### Applikation

Behandeling drinkwater

#### Water Safety Kit Microbiology Duo

In one hand or on the back: This test kit accompanies you in every situation with the complete microbiological test equipment. Test different sorts of bacteria directly on-site - with both incubators even simultaneously! Be sure of the quality of potable water everywhere, even when there is no laboratory nearby. The DI 20 incubator with its unique capability of heating and cooling as standard ensures reliable results.

#### User

Non Profit Organizations (NGO's) | Medical research and development | Institutes, Universities, Schools | World Health Organization (WHO)

#### Technische gegevens

**Afmetingen** 487 x 229 x 386 mm

**Gewicht met verpakking** 18.75 kg

## Leveromvang

- Marker Pen
- 2.5" Bailer with Rope
- Filter Flask, 500ml
- Tweezers
- Pump, 35ml stroke
- 5 Universal Container and spoon, 30 ml
- DB 513 Shaker Tube with lid, 100 ml
- 150 EZ Pack Filters
- 5 Pipette Disposable, 1.5 ml (EDT222)
- 2 DI 20 Incubator
- 200 Media Pads and dispenser
- Incubator Battery Lead/Acid
- Filtration Device
- 50 Petri Dish for DI 20 Incubator
- Lead Set IEC C13, (3 pin
- Membrane Lauryl Sulphate Broth, 40 g
- Instruction manual for Water Safety Kit Microbiology
- Instruction manual for DI 20 Incubator
- Case

## Accessoires

Omschrijving	Bestelnr.
Stroomkabelset	136300
Incubator Battery, Lead/Acid	192345
Flesjes	363590
Universal container with spoon and lid, 30 ml	366239
Plastic beker, gegradueerd, 100 ml, met deksel	385130
Petrischaaltjes, 10 stuks	400855
Rugzak voor Water Safety Kits	56A01400
Media Pads and Dispenser (Pack of 200)	56A016330
Maatbeker	56A018501
Membraanfiltratie	56A023801
Pincet	56A604201
Membraanfiltratie	56B000710
D004 TTC/MAC Dipslides (Macconkey)	56B010410
D005 TTC/Chromogeen Dipslides	56B010510
Membrane Lauryl Sulphate Broth, 38.1 gr.	56P069740
Pump (35 ml stroke)	56X610301
Pipette, disposable, 1.5 ml	BI-K10060
Marker Pen	BI-K14087

### Tintometer GmbH

Lovibond® Water Testing  
Schleefstraße 8-12  
44287 Dortmund  
Tel.: +49 (0)231/94510-0  
sales@lovibond.com  
www.lovibond.com  
Duitsland

### The Tintometer Limited

Lovibond House  
Sun Rise Way  
Amesbury, SP4 7GR  
Tel.: +44 (0)1980 664800  
Fax: +44 (0)1980 625412  
sales@lovibond.uk  
www.lovibond.com  
Verenigd Koninkrijk

### Tintometer China

9F, SOHO II C.  
No.9 Guanghualu,  
Chaoyang District,  
Beijing, 100020  
Customer Care China Tel.: 4009021628  
Tel.: +86 10 85251111 Ext. 330  
Fax: +86 10 85251001  
chinaoffice@tintometer.com  
www.lovibond.com  
China

### Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,  
Lebuh Nilam 2, Bandar Bukit Tinggi,  
Klang, 41200, Selangor D.E  
Tel.: +60 (0)3 3325 2285/6  
Fax: +60 (0)3 3325 2287  
lovibond.asia@tintometer.com  
www.lovibond.com  
Maleisië

### Tintometer Brazil

Caixa Postal: 271  
CEP: 13201-970  
Jundiaí – SP  
Tel.: +55 (11) 3230-6410  
sales@lovibond.us  
www.lovibond.com.br  
Brazilië

### Tintometer Inc.

6456 Parkland Drive  
Sarasota, FL 34243  
Tel: 941.756.6410  
Fax: 941.727.9654  
sales@lovibond.us  
www.lovibond.us  
Verenigde Staten

### Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2<sup>nd</sup>, 3<sup>rd</sup> & 4<sup>th</sup> Floor  
Sanathnagar Industrial Estate,  
Hyderabad, 500018  
Telangana  
Tel: +91 (0) 40 23883300  
Toll Free: 1 800 599 3891/ 3892  
indiaoffice@lovibond.in  
www.lovibondwater.in  
India

### Tintometer Spain

Postbox: 24047  
08080 Barcelona  
Tel.: +34 661 606 770  
sales@tintometer.es  
www.lovibond.com  
Spanje