

Lovibond® Water Testing Tintometer® Group



SD10 pH



- Transparent lid for easy maintenance and storage
- Ideal lid design for a stable stand
- Display colours & symbols for easy use
- New multi-parameter models for individual measurement requirements

Bestelnr.: 720010

Industrie

Gemeenten | Ngo | Openbaar zwembad | Overige industrieën | Privé-zwembad | Scheepvaart

Applikation

Afvalwaterzuivering | Behandeling drinkwater | Controle zwembadwater | Ketelwater | Koelwater

SD10 pH

With its waterproof housing, the Pocket Tester is ideal for regular checks in industrial plants or swimming pools.

Plus, the transparent lid not only provides a secure stand for the Pocket Tester in daily use, but also makes it easier to maintain the electrodes. A three-colour LED lighting system - blue, green and red - and new display icons make the tester easy to use and intuitive to operate. The use of energy saving LEDs and 4 x AAA batteries greatly extends the life of the new Pocket Tester range to up to 2000 hours.

With the included lanyard, each SD Series tester can be worn comfortably around the neck or attached to a pocket, waistband or belt with the clip on the head of the tester, allowing it to be taken anywhere without tiring the hands.

Meetbereik

Test Name	Meetbereik
pH	0 - 14
pH (mV)	-1000 - 1000
Temperature	0 - 50 °C

Technische gegevens

Nauwkeurigheid	pH: ± 0.01 (± 1 digit) pH (mV): ± 0.2 % FS Temperature: ± 0.5 °C
Display	3 coloured backlit LCD (blue: measurement, green: calibration, red: alarm)
Auto – OFF	Ja
IJking	1-, 2- of 3 puntsinstelling met automatische herkenning (NIST / IUPAC)
Stroomvoorziening	4 micro-batterijen (AAA)
Levensduur batterij	> 2000 h
Draagbaarheid	Draagbaar
Resolutie meetwaarde	pH: 0.01 ; pH (mV): 1
Resolution Temperature	0.1 °C
Calibration Points	1-3
Buffer Auto Recognition	USA: 1.68 / 4.00 / 7.00 / 10.01 / 12.45 pH NIST: 1.68 / 4.01 / 6.86 / 9.81 / 12.45 pH
Calibration Indicator	Ja
Stability Indicator	Ja
Auto-Hold Function	Ja
Alarm Limits	Ja
Compliance	CE
Afmetingen	40 x 40 x 180 mm
Gewicht met verpakking	268 g

Leveromvang

- Meter in a robust plastic case with hanger
- 4 batteries
- Lanyard
- pH 4, 7, 10 buffer tablets (1 Strip of 10 tablets each)
- 3 M KCl Storage Solution for pH/ORP
- Instruction manual
- Warranty declaration
- Certificate of Compliance

Omschrijving	Bestelnr.
pH-bufferoplossing 10,01 (25 °C) blauw, 90 ml, herleidbaar tot N.I.S.T.	721249
pH-bufferoplossing set pH 4, pH 7, pH 10 (25 °C), 90 ml	721250
pH-bufferoplossing 4,01 (25 °C) rood, 1 liter, herleidbaar tot N.I.S.T.	721252
pH-bufferoplossing 7,00 (25 °C), geel, 1 liter N.I.S.T.	721254
pH-bufferoplossing 10,01 (25 °C) blauw, 1 liter, herleidbaar tot N.I.S.T.	721256
3 M KCl-bewaringsoplossing voor pH/ORP, 25 ml (beschikbaar vanaf 2020)	726402
3 M KCl-bewaringsoplossing voor pH/ORP, 100 ml	726404

Accessoires

Omschrijving	Bestelnr.
Microbatterijen (AAA), set van 4	1950026
Maatbeker, 100 ml	384801
VE-water	461275
Buffertabletten pH 10	515600BT
Buffertabletten pH 10	515601BT
Buffertabletten pH 7	515610BT
Buffertabletten pH 7	515611BT
Buffertabletten pH 4	515620BT
Buffertabletten pH 4	515621BT
Vervangingselektrode voor SD10 pH	720015
pH-bufferoplossing 4,01 (25 °C): rood, 90 ml, herleidbaar tot N.I.S.T.	721247
pH-bufferoplossing 7,00 (25 °C) geel, 90 ml herleidbaar tot N.I.S.T.	721248

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Duitsland

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
Verenigd Koninkrijk

Tintometer China

9F, SOHO II C.
No.9 Guanghualu,
Chaoyang District,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 Ext. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
China

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Maleisië

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
Brazilië

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
Verenigde Staten

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
India

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
Spanje