



Lovibond® Water Testing

Tintometer® Group



Closed System Engineers Test Kit

The expert tool for reliable result



These Test Kits for cooling water monitoring are specially designed for engineers and mobile use. It provides the equipment and reagents required to monitor the correct amount and dosing levels of various important parameters. You can carry everything you need with you and test directly onsite. The Kit is provided in a robust carry case and includes all important equipment from measurement instruments to reagents.

- Complete analysis kit including molybdate, iron and conductivity
- Prevents corrosion and scale
- Accurate inhibitor dosage rates
- Reproducible, accurate results

Номер заказа: 56K701600

Промышленность

Другие отрасли | Коммуны | Нефтяная промышленность | Пищевая промышленность и производство напитков | Поставщики-энергии | Судоходство | Фармацевтическая промышленность | Химическая промышленность

Применение

Котельная вода | Обработка сточных вод | Охлаждающая вода

Closed System Engineers Test Kit

Due to its heat absorption capacity, water is the preferred cooling medium in many areas of industry. However, water always contains many dissolved and undissolved substances that can cause trouble even when travelling in closed cooling circuits. Extensive water treatment or chemical post-treatment by experts is always necessary. Problems are mainly caused by deposits, corrosion and microbiological activities. For the experts, a regular check of certain parameters is therefore just as essential as reliable equipment.

Диапазон измерений

Test Name	Диапазон измерений	Химический метод
Alkalinity M	50 - 2400 mg/L CaCO ₃	
Chloride	20 - 12000 mg/L Cl ⁻	
Copper HR T	0 - 5 mg/L Cu	Биквинолин
Iron LR T	0.05 - 1 mg/L Fe	PPST
Molybdate HR T	50 - 500 mg/L MoO ₄	Тиогликолят
Алюминий Т	0 - 0.3 mg/L Al	Эриохромоцианин Р
Гликоль	% PEG/MEG	
Жесткость общая	5 - 600 mg/L CaCO ₃	
Нитрит	10 - 2000 mg/L NaNO ₂	

Объем поставки

- Plastic case
- SD40 Multi Pocket Tester
- Checkit Comparator
- 50 - 500 mg/L Molybdate HR+ Checkit Disc
- 0 - 1.0 mg/L Iron LR Checkit Disc
- 0 - 0.3 mg/L Aluminium Checkit Disc
- 0 - 5.0 mg/L Copper HR Free and Total Checkit Disc
- 2 Plastic vials 10 mL with lid
- Aluminium No.1 tablets
- Aluminium No.2 tablets
- Copper No. 1 Test tablets
- Copper No. 2 Test tablets
- Molybdate No.1 tablets
- Molybdate No.2 tablets
- pH Buffer 4 solution
- pH Buffer 7 solution
- pH Buffer 10 solution
- Iron Reagent FE5
- Hardness Total Buffer TH2
- Hardness Total Indicator TH1P
- Hardness LR Titrant TH3
- Hardness HR Titrant TH4
- 1413µS Standard Conductivity Solution
- Alkalinity 4.5 Indicator TA4
- Alkalinity LR Titrant TA3
- Alkalinity HR Titrant PA2/TA2
- Nitrite Indicator N1
- Nitrite Titrant N2
- Chloride Indicator CC1
- Chloride LR Titrant CC2
- Chloride HR Titrant BC2
- 2 Titration jars with caps, plastic, 60 mL
- Syringe plastic, 20 mL
- Wash Bottle, 250 mL
- Syringe Filter Holder 25 mm
- Filter Circle 0.45 µm, 25 mm
- Glycol Refractometer, RHA-218ATC
- Instruction manual

Аксессуары

Заголовок	Номер заказа
Компаратор СНЕСКИТ®	145000
Шприц, пластик, 20 мл	56A006501
Банка для титрования с крышкой, пластик, 60 мл	56A006701
Пластиковый шприц, 5 мл	56A008501
Промывная склянка 250 мл	56A009001
Пластиковый шприц, 1 мл	56A013501
Фильтровальный диск 0,45 мкм, 25 мм	56A020050
Пробирка 5/10 мл + CAP	56A600401
Рефрактометр для гликоля	56M000401

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Германия

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
ВЕЛИКОБРИТАНИЯ

Tintometer China

Room 1001, China Life Tower
16 Chaoyangmenwai Avenue,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 App. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
Китай

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Малайзия

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
Бразилия

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
США

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
Индия

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
Испания

Возможны технические изменения

Отпечатано в Германии

Lovibond® and Tintometer® are Trademarks of the Tintometer Group of Companies