

# Lovibond® Water Testing

## Tintometer® Group



## Pooltester Klor LR/pH karton içinde

The fast assistants for Pool Testing



- Easy to use
- RAPID tablets fast dissolving
- Highest accuracy

Ürün No: 151613

### High measuring accuracy

Results are shown with high measurement accuracy. They have a certainty and can counteract if necessary.

### Intuitive

Self-explanatory is the operation of the tester. And if you are ever unsure: A look at the packaging and the intuitively understandable instructions are sufficient to eliminate the ambiguity.

### Ergo mixed

The ergonomically optimized lid ensures a secure grip.

### Sanayi

Yerel Havuz

### Uygulama

Afet Yönetimi | Dezenfeksiyon Kontrol | Havuz Suyu Kontrol | İçme-Suyu Arıtma | Soğutma Suları | Yiyecek ve İçecek

### Pooltester Klor LR/pH karton içinde

Sürdürülebilir bir geleceğe doğru bir adım: Pooltester Cl LR, pH artık geri dönüştürülebilir ambalajlarda mevcuttur.

Pooltester klor içeriği ve pH değerinin eşzamanlı olarak ölçülmesini sağlar: havuz suyunu doldurun, Rapid tabletleri ekleyin ve renk değişimini ölçekle karşılaştırın. Bu, havuz sahiplerinin havuzlarının su bakımını takip etmelerini kolaylaştırır.

### Ölçüm Aralığı

| Test Name | Ölçüm Aralığı                  |
|-----------|--------------------------------|
| Chlorine  | 0.1 - 3.0 mg/L Cl <sub>2</sub> |
| pH        | 6.8 - 8.2 pH                   |

# Teslimat Kapsam Metni

- Pooltester
- Tablet reagents for 20 tests
- Instruction manual

## **Tintometer GmbH**

Lovibond® Water Testing  
Schleefstraße 8-12  
44287 Dortmund  
Tel.: +49 (0)231/94510-0  
sales@lovibond.com  
www.lovibond.com  
Almanya

## **The Tintometer Limited**

Lovibond House  
Sun Rise Way  
Amesbury, SP4 7GR  
Tel.: +44 (0)1980 664800  
Fax: +44 (0)1980 625412  
sales@lovibond.uk  
www.lovibond.com  
Birleşik Krallık

## **Tintometer China**

Room 1001, China Life Tower  
16 Chaoyangmenwai Avenue,  
Beijing, 100020  
Customer Care China Tel.: 4009021628  
Tel.: +86 10 85251111 App. 330  
Fax: +86 10 85251001  
chinaoffice@tintometer.com  
www.lovibond.com  
Çin

## **Tintometer South East Asia**

Unit B-3-12, BBT One Boulevard,  
Lebuh Nilam 2, Bandar Bukit Tinggi,  
Klang, 41200, Selangor D.E  
Tel.: +60 (0)3 3325 2285/6  
Fax: +60 (0)3 3325 2287  
lovibond.asia@tintometer.com  
www.lovibond.com  
Malezya

## **Tintometer Brazil**

Caixa Postal: 271  
CEP: 13201-970  
Jundiaí – SP  
Tel.: +55 (11) 3230-6410  
sales@lovibond.us  
www.lovibond.com.br  
Brezilya

## **Tintometer Inc.**

6456 Parkland Drive  
Sarasota, FL 34243  
Tel: 941.756.6410  
Fax: 941.727.9654  
sales@lovibond.us  
www.lovibond.us  
Amerika Birleşik Devletleri

## **Tintometer India Pvt. Ltd.**

Door No: 7-2-C-14, 2<sup>nd</sup>, 3<sup>rd</sup> & 4<sup>th</sup> Floor  
Sanathnagar Industrial Estate,  
Hyderabad, 500018  
Telangana  
Tel: +91 (0) 40 23883300  
Toll Free: 1 800 599 3891/ 3892  
indiaoffice@lovibond.in  
www.lovibondwater.in  
Hindistan

## **Tintometer Spain**

Postbox: 24047  
08080 Barcelona  
Tel.: +34 661 606 770  
sales@tintometer.es  
www.lovibond.com  
İspanya