

# Lovibond® Water Testing Tintometer® Group



## Water Safety Kit Combined Duo

Extended special test kit for potable water not only for NGOs



- Emergency response kits for disaster situations and water monitoring outback
- Equipment for chemical and microbiological analysis in two cases

Ürün No: 56K681252

You need safe, fast and reliable test results in difficult situations like disasters or potable water monitoring in rural environments without laboratories? You immediate certainty about the quality of potable water to save health and lives? Take all you need with you and combine microbiological methods and simple chemical methods of analysis. Take the smaller Water Safety Kit Chemical in one hand and the bigger Water Safety Kit Microbiology! on your back in the optional backpack. You are able to test different sorts of bacteria on-site due to the heating and cooling function of our DI 20 petri dish incubator.

### Sanayi

Belediyeler | Deniz Sanayi | Diğer Sanayiler | Güç ve Enerji | Halka-Açık Havuz | Havuz Hizmetleri | İlaç Sanayi | Kimya Sanayi | Sivil-Toplum Kuruluşları | Yağ ve Gaz | Yerel Havuz | Yiyecek ve İçecek-Sanayi

### Uygulama

Atık Su Arıtma | Dezenfeksiyon Kontrol | Ham Su Arıtma | Havuz Suyu Kontrol | İçme Suyu Arıtma | Kazan Suları | Soğutma Suları

### Water Safety Kit Combined Duo

Combine both microbiological and physio chemical analysis of some of the most critical parameters to determine the suitability of water for drinking. Take the extended complete equipment including 2 incubators for bacteria tests with you everywhere you need immediate and reliable test result.

With their rugged, waterproof cases the kits are easy to transport. Even non technical personal is able to use it due to pictorial instructions.

### User

Non Profit Organizations (NGO's) | Medical research and development | Institutes, Universities, Schools | World Health Organization (WHO)

## Ölçüm Aralığı

Test Name	Ölçüm Aralığı	Kimyasal Metod
Aerobik bakteri (TVC)/ E.Coli, koliformlar		
Aerobik bakteriler (TVC)/ koliformlar		
Amonyum T	0.02 - 1 mg/L N	İndofenol Mavisi
Chlorine T	0.01 - 6.0 mg/L Cl <sub>2</sub> <sup>a)</sup>	DPD
Conductivity	0 - 20 mS/cm	
Nitrat T	0.08 - 1 mg/L N	Çinko İndirgemesi / NED
pH	0 - 14	
Turbidity	0.01 - 1100 NTU	

## Teslimat Kapsam Metni

- 1 case Water Safety Kit Microbiology Duo
- 1 case Water Safety Kit Chemical

## Aksesuarlar

Başlık	Ürün No
Güç kablosu seti	136300
Incubator Battery, Lead/Acid	192345
Satz Trübungsstandards T-CAL (<0,1, 20, 200, 800 NTU)	194150
Şişeler	363590
Universal container with spoon and lid, 30 ml	366239
Plastik kap, 100 ml derecelendirmeli, kapaklı	385130
Plastik kap, 30 ml derecelendirmeli, kapaklı	385132
Petri kapları, 10 adet.	400855
ValidCheck Klor 1,5 mg/l	48105510
Karıştırma çubuğu ve toz kaşığı	56A006601
Su Güvenlik Kitleri için sırt çantası	56A01400
Media Pads and Dispenser (Pack of 200)	56A016330
Ölçü kabı	56A018501
Diyaframlı filtreleme	56A023801
Cımbız	56A604201
Diyaframlı filtreleme	56B000710
D004 TTC/MAC daldırma kızakları (Macconkey)	56B010410
D005 TTC/Kromojenik daldırma kızakları	56B010510
KS3-pH 10 tampon çözelti 65 ml	56L000365
KS6-pH 4 tampon çözelti 65 ml	56L000665
KS9-pH 7 tampon çözelti 65 ml	56L000965
KS16-1413US standart iletkenlik çözeltisi	56L001665
Membrane Lauryl Sulphate Broth, 38.1 gr.	56P069740
Pump (35 ml stroke)	56X610301
Pipette, disposable, 1.5 ml	BI-K10060
Marker Pen	BI-K14087

### Tintometer GmbH

Lovibond® Water Testing  
Schleefstraße 8-12  
44287 Dortmund  
Tel.: +49 (0)231/94510-0  
sales@lovibond.com  
www.lovibond.com  
Almanya

### The Tintometer Limited

Lovibond House  
Sun Rise Way  
Amesbury, SP4 7GR  
Tel.: +44 (0)1980 664800  
Fax: +44 (0)1980 625412  
sales@lovibond.uk  
www.lovibond.com  
Birleşik Krallık

### Tintometer China

Room 1001, China Life Tower  
16 Chaoyangmenwai Avenue,  
Beijing, 100020  
Customer Care China Tel.: 4009021628  
Tel.: +86 10 85251111 App. 330  
Fax: +86 10 85251001  
chinaoffice@tintometer.com  
www.lovibond.com  
Çin

### Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,  
Lebuh Nilam 2, Bandar Bukit Tinggi,  
Klang, 41200, Selangor D.E  
Tel.: +60 (0)3 3325 2285/6  
Fax: +60 (0)3 3325 2287  
lovibond.asia@tintometer.com  
www.lovibond.com  
Malezya

### Tintometer Brazil

Caixa Postal: 271  
CEP: 13201-970  
Jundiaí – SP  
Tel.: +55 (11) 3230-6410  
sales@lovibond.us  
www.lovibond.com.br  
Brezilya

### Tintometer Inc.

6456 Parkland Drive  
Sarasota, FL 34243  
Tel: 941.756.6410  
Fax: 941.727.9654  
sales@lovibond.us  
www.lovibond.us  
Amerika Birleşik Devletleri

### Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2<sup>nd</sup>, 3<sup>rd</sup> & 4<sup>th</sup> Floor  
Sanathnagar Industrial Estate,  
Hyderabad, 500018  
Telangana  
Tel: +91 (0) 40 23883300  
Toll Free: 1 800 599 3891/ 3892  
indiaoffice@lovibond.in  
www.lovibondwater.in  
Hindistan

### Tintometer Spain

Postbox: 24047  
08080 Barcelona  
Tel.: +34 661 606 770  
sales@tintometer.es  
www.lovibond.com  
İspanya