

# Lovibond® Water Testing Tintometer® Group



## Water Safety Kit Chemical

Mobile Laboratory for chemical parameters of potable water testing



- Analysis of most important chemical parameters
- Complete test equipment in conjunction with Water Safety Kits Microbiology
- No further test equipment required for chemical analysis

Ürün No: 56K681253

You need safe, fast and reliable test results in difficult situations like disasters or potable water monitoring in rural environments without laboratories? The Lovibond® Water Safety Kit Chemical offers the most important methods and equipment in one case, designed to accompany you in each situation.

It comprises our MD 600 photometer, TB 211 IR turbidity meters and SD 50 AND SD 70 pocket electrochemistry meters. This kit can be used in conjunction with the Water Safety Kits Microbiology to provide a complete suite of tests for potability of water.

The kit comes with reagents for Chlorine, Ammonia, Nitrite, Nitrate, pH, Conductivity, Turbidity but can be expanded using our reagent packs to include any of the tests for the MD 600.

### Sanayi

Belediyeler | Deniz Sanayi | Diğer Sanayiler | Güç ve Enerji | Halka-Açık Havuz | Havuz Hizmetleri | İlaç Sanayi | Kimya Sanayi | Sivil-Toplum Kuruluşları | Yağ ve Gaz | Yerel Havuz | Yiyecek ve İçecek-Sanayi

### Uygulama

Atık Su Arıtma | Dezenfeksiyon Kontrol | Ham Su Arıtma | Havuz Suyu Kontrol | İçme Suyu Arıtma | Kazan Suları | Soğutma Suları

### Water Safety Kit Chemical

Take everything you need for chemical analysis of potable water quality in one hand. This Water Safety Kit is the smallest in the range, portable in a rugged, waterproof case with 200 reagents for each parameter, allowing numerous measurements for its small footprint.

It is easy to use even by nonprofessional operators due to simple pictorial instructions.

### User

Non Profit Organizations (NGO's) | Medical research and development | Institutes, Universities, Schools | World Health Organization (WHO)

## Ölçüm Aralığı

| Test Name    | Ölçüm Aralığı                                 | Kimyasal Metod          |
|--------------|---|-------------------------|
| Amonyum T    | 0.02 - 1 mg/L N                               | Indofenol Mavisi        |
| Conductivity | 0 - 20 mS/cm                                  |                         |
| Klor T       | 0.01 - 6.0 mg/L Cl <sub>2</sub> <sup>a)</sup> | DPD                     |
| Nitrat T     | 0.08 - 1 mg/L N                               | Çinko İndirgemesi / NED |
| pH           | 0 - 14  |                         |
| Turbidity    | 0.01 - 1100 NTU                               |                         |

## Teknik Veri

|                   |                    |
|-------------------|--------------------|
| Boyutlar          | 411 x 168 x 322 mm |
| Ambalajlı Ağırlık | 4.5 kg             |

## Aksesuarlar

| Başlık  | Ürün No   |
|---|-----------|
| Satz Trübungsstandards T-CAL (<0,1, 20, 200, 800 NTU) | 194150    |
| Plastik kap, 100 ml derecelendirmeli, kapaklı         | 385130    |
| Plastik kap, 30 ml derecelendirmeli, kapaklı          | 385132    |
| ValidCheck Klor 1,5 mg/l                              | 48105510  |
| Karıştırma çubuğu ve toz kaşığı                       | 56A006601 |
| Su Güvenlik Kitleri için sırt çantası                 | 56A01400  |
| KS3-pH 10 tampon çözelti 65 ml                        | 56L000365 |
| KS6-pH 4 tampon çözelti 65 ml                         | 56L000665 |
| KS9-pH 7 tampon çözelti 65 ml                         | 56L000965 |
| KS16-1413US standart iletkenlik çözeltisi             | 56L001665 |

## Teslimat Kapsam Metni

- MD 600 Photometer
- TB 211 IR Turbidimeter
- SD 70 Conductivity Meter
- SD 50 pH Meter
- KS16 (65 ml)
- KS9 (65 ml)
- KS6 (65 ml)
- KS3 (65 ml)
- 250 DPD No.1 Photometer tablet reagents
- 250 DPD No.3 Evo Photometer tablet reagents
- Nitrate Test Powder R282, 15 g
- 100 Nitrate Test tablet reagents
- Nitrate Test Tube Cali., 20 ml
- 250 Nitrite LR tablet reagents
- 250 Ammonia No.1 tablet reagents
- 250 Ammonia No.2 tablet reagents
- Shaker Tube Plastic, 30 ml
- DB513 Shaker Tube with lid, 100 ml
- Stirring Rod & Powder Scoop
- 3 Round Cuvette 24 mm with lid and seal
- 2 Round Cuvette 24 mm x 54 mm
- 2 Lid 24 mm, black
- T-CAL Turbidity Standard < 0.1 NTU
- T-CAL Turbidity Standard 20 NTU
- T-CAL Turbidity Standard 200 NTU
- T-CAL Turbidity Standard 800 NTU
- Instruction manual for Water Safety Kit Chemical
- Case

---

**Tintometer GmbH**

Lovibond® Water Testing  
Schleefstraße 8-12  
44287 Dortmund  
Tel.: +49 (0)231/94510-0  
sales@lovibond.com  
www.lovibond.com  
Almanya

**The Tintometer Limited**

Lovibond House  
Sun Rise Way  
Amesbury, SP4 7GR  
Tel.: +44 (0)1980 664800  
Fax: +44 (0)1980 625412  
sales@lovibond.uk  
www.lovibond.com  
Birleşik Krallık

**Tintometer China**

Room 1001, China Life Tower  
16 Chaoyangmenwai Avenue,  
Beijing, 100020  
Customer Care China Tel.: 4009021628  
Tel.: +86 10 85251111 App. 330  
Fax: +86 10 85251001  
chinaoffice@tintometer.com  
www.lovibond.com  
Çin

**Tintometer South East Asia**

Unit B-3-12, BBT One Boulevard,  
Lebuh Nilam 2, Bandar Bukit Tinggi,  
Klang, 41200, Selangor D.E  
Tel.: +60 (0)3 3325 2285/6  
Fax: +60 (0)3 3325 2287  
lovibond.asia@tintometer.com  
www.lovibond.com  
Malezya

**Tintometer Brazil**

Caixa Postal: 271  
CEP: 13201-970  
Jundiaí – SP  
Tel.: +55 (11) 3230-6410  
sales@lovibond.us  
www.lovibond.com.br  
Brezilya

**Tintometer Inc.**

6456 Parkland Drive  
Sarasota, FL 34243  
Tel: 941.756.6410  
Fax: 941.727.9654  
sales@lovibond.us  
www.lovibond.us  
Amerika Birleşik Devletleri

**Tintometer India Pvt. Ltd.**

Door No: 7-2-C-14, 2<sup>nd</sup>, 3<sup>rd</sup> & 4<sup>th</sup> Floor  
Sanathnagar Industrial Estate,  
Hyderabad, 500018  
Telangana  
Tel: +91 (0) 40 23883300  
Toll Free: 1 800 599 3891/ 3892  
indiaoffice@lovibond.in  
www.lovibondwater.in  
Hindistan

**Tintometer Spain**

Postbox: 24047  
08080 Barcelona  
Tel.: +34 661 606 770  
sales@tintometer.es  
www.lovibond.com  
İspanya